


Second Harmonic Generation of the FCPA μ Jewel

This note describes a simple set-up for second-harmonic generation (also known as “frequency-doubling”) of an FCPA μ Jewel to obtain green, pulsed light using off-the-shelf components.



522 nm beam profile obtained by frequency doubling an IMRA FCPA μ Jewel laser




Figure 1 LBO crystal SHG conversion with singlet lens

Set-up

Focus the laser output beam into the LBO crystal (See Figure 1). It is not critical that the focus be perfectly centered in the crystal. Recommended optics are given in the next section. The beam should be incident nearly perpendicular to the crystal front surface. However, a minor tilt adjustment may be necessary to get the best conversion efficiency. While there are other suppliers for LBO mounted in an oven for temperature control, a product from Newlight Photonics is described in the table below for your convenience. See Table 1 and Figure 2 for more information.

LBO temperature setting: 174.4°C (SNLO- see reference), adjusted for the best conversion

Among the various nonlinear crystals available (BBO, KTP, BiBO, KNbO₃, PPLN, etc.), we recommend the use of Lithium Triborate (LBO) for the following reasons:

- Non-Critical Phase Matching (NCPM) condition
 - No walk-off
 - Large acceptance angle → ease of alignment

Since the energy per pulse from the FCPA μ Jewel changes for different repetition rates, different focusing conditions are required for optimal conversion.

- Small Group Velocity Mismatch (GVM) allows longer crystal to be used → greater SHG output
- High damage threshold
Intensity = 45 GW/cm², t = 100 ps @ I = 532 nm (ref. Newlight Photonics)

LBO Crystal 3 mm thick, 6 mm square face Cut angles: $\theta = 90^\circ$ $\phi = 0^\circ$ 1045 nm to 522 nm conversion (NCPM) AR coating: 1045 and 522 nm	Newlight Photonics
Oven & Precision Temperature Controller	Newlight Photonics
2" Kinematic Mirror Mount (for oven)	Thorlabs KM 200

Table 1 Crystal assembly components
www.newlightphotonics.com www.thorlabs.com



Figure 2 Equipment Photos

IMRA America, Inc.
1044 Woodridge Ave.
Ann Arbor, MI 48105

Main: (734) 930-2560
Fax: (734) 930-9957
lasers@imra.com
www.imra.com

IMRA The Femtosecond Fiber Laser Company

Optics Selection

The focusing lens should be anti-reflection coated for 1 μm light and the collimating lens should be anti-reflection coated for $\sim 520\text{ nm}$ to obtain optimal results

Focusing Lens: The FCPA μJewel pulse energy varies for different repetition rates, so different focusing optics will produce the optimum conversion efficiency. The lenses in Table 2 are based on typical values of pulse energy at different repetition rates.

repetition rate	focusing lens
100 kHz	f = 400-500 mm
200 kHz	f = 300-500 mm
500 kHz	f = 200-300 mm
1 MHz	f = 150-250 mm
5 MHz	f = 100 mm

Table 2 Recommended focusing lenses for different repetition rates

Collimating lens: Selected for the desired beam size

Conversion Efficiency

The conversion efficiency peaks at a certain crystal temperature (see Figure 3). The conversion efficiency can be up to 50%, as shown in Figure 4. It is recommended that maximum input laser power is used to generate SHG. If less power at $\sim 520\text{ nm}$ is required, attenuate the converted beam.




Figure 3 SHG power vs. temperature efficiency at 200 kHz. Optimum is at 174.8 °C




Figure 4 SHG efficiency at 200 kHz and f = 200 mm

Spectra

The spectra for the input fundamental and the generated second harmonic are shown below.




Figure 5 FCPA μJewel D-400 fundamental (top) and second harmonic (bottom) spectra

Summary

Using off-the-shelf optics and a simple set-up, conversion of the 1 μm output of the FCPA μJewel laser to green output of $\sim 520\text{ nm}$ is achieved with conversion efficiency of up to 50%.