

## IMRA Product Information

### Femtolite HX-150 Wavelength Selection Settings

<b>Filter Wheel Settings</b>	<b>F</b>	<b>F</b>	<b>F</b>	<b>H</b>	<b>H</b>	<b>H</b>	<b>F + H</b>	<b>F + H</b>	<b>F + H</b>
<b>Slider Setting</b>	F	H	Between <b>F</b> and <b>H</b>	F	H	Between <b>F</b> and <b>H</b>	F	H	Between <b>F</b> and <b>H</b>
<b>1620 nm Band</b>	—	—	—	Min	Max	Adjustable	Min	Max	Adjustable
<b>810 nm Band</b>	Max	Min	Adjustable	—	—	—	Max	Min	Adjustable

